



1-877-699-8439

Finger Pulse Oximeter from Oxy-View



Quick and easy
spot check of heart rate
and oxygen saturation

Six (6) colorful and easy
to read display options

Two (2) different
waveform displays

Display values especially
clear and easy to read in
low ambient light situations

More than 30 hours of
continuous operation

Automatic power off

Low battery indicator

Uses 2-AAA batteries

Full 2-year
limited warranty

FDA approved

CE 0123
IEC 61601-1-1

Product and price subject to change to reflect
continuing product improvements, changes in
specifications and availability.

Finger Pulse Oximeter



Your vital signs: Temperature...Pulse...Respirations...Blood Pressure

For over 140 years there have been just four measurements we call "Vital Signs." Finally, due to the combination of 21st century technology and healthy competition, the 5th vital sign...finger pulse oximetry is available at price never seen before. **Why pay more?**

Order Number: OX300C2H

- Display: Color OLED
Parameters: SpO₂, PR, Pleth bar, Pleth waveform
Information: Low battery indicator
Direction of the display can be changed
- SpO₂: Range: 35-99%
Accuracy: 70% -100% ± 2%
0% - 69% unspecified
Method: Dual wavelength LED
- Pulse Rate: Range: 30 - 240 BPM
Accuracy: ± 2 BPM or 2%
Resolution: 1 BPM
- Battery: Type: 2-AAA Alkaline or rechargeable batteries
- Environmental: Operating Temperature: 5° C to 40° C (41°F to 104°F)
Storage Temperature: -10° C to 40° C (14° to 104° F)
Operating Humidity: 15%RH to 95% RH
Storage Humidity: 10%RH to 95% RH



109 Inverness Drive East, Suite C
Englewood, Colorado 80112-5105

www.oxyview.com Phone: 877-699-8439 Fax: 303-790-4588